## ASPINA - DO M

- EN User manual
- DE Benutzerhandbuch
- FR Manuel de l'utilisateur
- RU Руководство пользователя
- PL Instrukcja obsługi
- SK Návod na použitie
- CS Návod k použití

CE 2460

THE DENTAL ASPIRATOR  
DENTALE ABSAUGEINHEIT  
L'ASPIRATEUR DENTAIRE  
ДЕНТАЛЬНОЕ ОТСАСЫВАЮЩЕЕ УСТРОЙСТВО  
ODSYSARKA DENTALNA  
DENTÁLNA ODSÁVAČKA  
DENTÁLNÍ ODSÁVAČKA

ASPINA - DO M

**EKOM spol. s r. o.**  
Priemyselná 5031/18  
SK-921 01 Piešťany  
Slovak Republic  
tel.: +421 33 7967255  
fax: +421 33 7967223

[www.ekom.sk](http://www.ekom.sk)  
email: [ekom@ekom.sk](mailto:ekom@ekom.sk)

DATE OF LAST REVISION  
DATUM DER LETZTEN ÜBERARBEITUNG  
DATE DE MISE A JOUR  
ДАТА ПОСЛЕДНЕГО ПЕРЕСМОТРА  
DATA WYDANIA  
DÁTUM POSLEDNEJ REVÍZIE  
DATUM POSLEDNÍ REVIZE

10/2020

NP-DO M-37\_10-2020  
112000080-000

<b>EN</b>	<b><u>CONTENTS..... 5</u></b>
<b>DE</b>	<b><u>INHALT ..... 22</u></b>
<b>FR</b>	<b><u>SOMMAIRE..... 40</u></b>
<b>RU</b>	<b><u>ОГЛАВЛЕНИЕ..... 58</u></b>
<b>PL</b>	<b><u>SPIS TREŚCI ..... 76</u></b>
<b>SK</b>	<b><u>OBSAH ..... 93</u></b>
<b>CS</b>	<b><u>OBSAH .....110</u></b>

**CONTENTS**

<b>IMPORTANT INFORMATION</b> .....	<b>6</b>
1. WARNINGS.....	6
1.1. CE Marking.....	6
1.2. General warnings.....	6
1.3. General safety warnings.....	6
1.4. Safety warnings regarding the protection against electric current.....	6
1.5. Alert notices and symbols.....	7
2. PRODUCT INFORMATION.....	8
2.1. Intended Use.....	8
2.2. Contraindications and side-effects.....	8
2.3. Product Description.....	8
3. TECHNICAL DATA.....	9
3.1. FAD efficiency correction for differences in elevation.....	9
3.2. Electromagnetic compatibility declaration.....	10
4. FUNCTION DESCRIPTION.....	13
4.1. Function Description.....	13
4.2. Detailed Description of the Suction Part Function.....	13
<b>ASSEMBLY</b> .....	<b>13</b>
5. STORAGE AND TRANSPORT CONDITIONS.....	13
6. INSTALLATION AND INITIAL START-UP.....	14
6.1. Environmental Conditions.....	14
6.2. Electric Connection.....	14
6.3. Initial Start-up.....	14
<b>USE</b> .....	<b>15</b>
7. OPERATION.....	15
7.1. Switching on the Dental Aspirator.....	15
7.2. Use of antifoaming tablets.....	15
7.3. Filter use in the separation vessel.....	15
8. MAINTENANCE INTERVALS - USER / TECHNICIAN.....	16
9. MAINTENANCE, CLEANING AND DISINFECTION.....	16
9.1. Inlet Sieve Cleaning.....	17
9.2. Lubricating seals and sliding closures.....	17
9.3. Disinfection of Tubings and Separation Vessel.....	17
9.4. Cleaning and disinfection of the exterior surfaces of the product.....	17
9.5. Replacement of the Output Pre-filter.....	17
9.6. Replacement of Output Filter.....	18
9.7. Cannulas handling and preparation operation.....	18
9.8. Replacing the filter in the separation vessel.....	19
10. PUTTING OUT OF OPERATION.....	19
<b>DISPOSAL</b> .....	<b>19</b>
11. DISPOSAL OF APPLIANCE.....	19
<b>TROUBLESHOOTING</b> .....	<b>20</b>
12. INFORMATION ON REPAIR SERVICE.....	20
13. SOLVING COMMON PROBLEMS.....	20
<b>ANNEXES</b> .....	<b>126</b>
14. WIRING DIAGRAMS.....	126
15. FUNCTION DIAGRAM.....	128
16. ENCLOSURE NO. 1.....	129
17. INSTALLATION RECORD.....	130

## **IMPORTANT INFORMATION**

### **1. WARNINGS**

#### **1.1. CE Marking**

Products marked with **CE** mark of compliance meet requirements of European Union directives (93/42/EEC).

#### **1.2. General warnings**

- The installation, operation and maintenance manual is an integral part of the appliance. It is necessary to always keep this document close to the appliance. Strict observance of this manual is a prerequisite for the correct operation of the appliance.  
The safety of operating personnel and failure-free operation of the appliance are ensured only when using the original components of the appliance. It is possible to use accessories and spare parts mentioned in the technical documentation or clearly permitted by the manufacturer.  
When used with non authorized accessories or consumable material, the manufacturer cannot assume responsibility for the safe operation and functionality of the device.
- The Guarantee does not cover damages that originate due to the use of non authorized accessories or consumable material other than those recommended by the manufacturer.
- The manufacturer assumes responsibility regarding safety, reliability and function only if
  - The installation, calibration, amendments, extensions and repairs are made by the manufacturer or his representative or a service organization authorized by the manufacturer,
  - The appliance is used in accordance with the installation, operation and maintenance manual.
- The installation, operation and maintenance manual, at the time of printing, corresponds to the design of the appliance and its state according to the relevant safety and technical standards. The manufacturer reserves all copyrights for the given wiring schemes, methods and names.
- This user manual is the original instructions. Translation is performed in accordance with the best available knowledge.

#### **1.3. General safety warnings**

The manufacturer developed and constructed the appliance so that damage would not occur when the appliance is used for its intended purpose. The manufacturer considers it his obligation to describe the following safety measures in order to avoid further damages.

- When operating the appliance, it is necessary to observe laws and regional regulations valid in the place of usage. In order to ensure safe course of works, the operator and user are responsible for the observation of regulations.
- The original packaging should be kept for the possible return of the unit. Only original packaging guarantees an optimal protection of the appliance during transportation. If it would be necessary to return the appliance during warranty period, the manufacturer is not responsible for damages caused by incorrect packing.
- It is necessary that the user ensures the appliance is safe to use prior to usage.
- The user must familiarize himself with the correct operation of appliance.
- If an undesirable event occurs in the operation of appliance, the user is obliged to immediately inform his supplier to this event.
- This product is not intended for use in areas with the risk of explosion.


#### **1.4. Safety warnings regarding the protection against electric current**

- The appliance must only be connected to an appropriate power source that has correct grounding.
- Prior the connecting the compressor, verify whether the mains voltage and frequency specified on the equipment are in accordance with the local supply.
- Prior to putting into operation, check for possible damages on the appliance and the air connectors. Damaged cables and sockets/plugs must be replaced immediately.

- In the case of a dangerous situation or a technical failure, immediately disconnect the appliance from mains supply.
- During all repairs and maintenance:
  - ensure that the mains plug is removed from the power socket
  - pressure pipes must be air vented
  - pressure must be released from pressure tank.
- This appliance can only be installed only by a qualified expert.

**1.5. Alert notices and symbols**

For your information, the symbols below are used in the installation, operation and maintenance manual, on packaging materials and on the product :

	General warnings
	Warning
	Danger, electric shock hazard.
	Read the operating instructions
	Refer to instruction manual
	CE – marking
	Attention! Hot surface.
	Handling mark on package – FRAGILE
	Handling mark on package – THIS SIDE UP
	Handling mark on package – KEEP DRY
	Handling mark on package – TEMPERATURE LIMITATIONS
	Handling mark on package – LIMITED STACKING
	Mark on package – RECYCLABLE MATERIAL
	Connection of protective ground wire
	Fuse
	Danger of biological hazard
	Indicates a medical device that has not undergone sterilization

Sterilizable in a steam sterilizer (autoclave) at specified temperature

Manufacturer

## 2. PRODUCT INFORMATION

### 2.1. Intended Use

The dental aspirator ASPINA DO M is a mobile equipment extending the possibilities of procedures provided for by a physician and ensuring the increase of ergonomics at his work. It is suitable for the dental units not equipped with the suction equipment and a separator. It is intended for the dental practice needs for aspiration, separation and entrapping of the waste into the built-in separation vessel. With regard to its mobility it enables a simple transport between the workplaces and can be used above all where the connection to a sewerage system is impossible.

The dental aspirator is designed for the operation in dry, ventilated rooms, with the ambient temperature ranging between +5°C up to +40°C, and relative humidity not exceeding the value of 70%.

- The dental aspirator may not be exposed to rain. The equipment may not be operated in the moist or wet environs. In addition, avoid the use near gases or combustible liquids.
- Other use, or the use exceeding this scope, can not be regarded as the intended use. The manufacturer is not responsible for the resulting damages. The risk will be borne exclusively by the operator / user.

### 2.2. Contraindications and side-effects

There are no contraindications or side-effects known.

### 2.3. Product Description

(Annex No.1)

The mobile dental aspirator ASPINA DO M is built-up on a movable bogie on which is situated a box, lined by a noise damping material. Inside the box there is situated a fan cooled suction aggregate (9) with the electric distribution and the waste entrapping separation vessel (11). In the lower part – under the box – there is situated a silencer with an output filter (14) and a pre-filter (15) ensuring the air filtration from the suction aggregate. In the upper – narrowed part of the aspirator - there is situated a holder for suction tubings (2) equipped with cannulas (1), separation automatics and terminal box with fuses. On the lateral part there is situated a main switch (5), over which are placed the indicators for the network (3) and for the state of filling of the separation vessel (4).

**3. TECHNICAL DATA**

		<b>DO M</b>
Rated voltage / frequency	V / Hz	230 / 50; 230 / 60*; 110 / 60*
Nominal current	A	2.6
maximum flowrate :	l/min	1100
maximum underpressure	kPa	12
Sound level	dB(A)	≤48
Mode of operation of the aspirator		permanent S 1
Dimensions of the aspirator	mm	565x350x860
Dimensions of the aspirator with carton (h x w x d)	mm	950x560x595
Weight of the aspirator	kg	36
Weight of the aspirator with Carton	kg	40
Model according to STN EN 60 601-1 (IEC 60601-1)		equipment type B, class I
Classification pursuant to MDD 93/42/EEC, 2007/47/EC		Ila

(\*) The rated voltage and frequency can be change by special order

**Climatic conditions for the operation:**

**Temperature** +5°C up to 40°C

**Relative air humidity** up to +70%

**3.1. FAD efficiency correction for differences in elevation**

*FAD correction table*

Elevation [mamsl]	0 - 1500	1501 - 2500	2501 - 3500	3501 - 4500
FAD [l/min]	FAD x 1	FAD x 0.8	FAD x 0.71	FAD x 0.60

FAD efficiency refers to conditions at an elevation of 0 mamsl:

Temperature: 20°C

Atmospheric pressure: 101325 Pa

Relative humidity: 0%

### 3.2. Electromagnetic compatibility declaration

Medical equipment needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided below.

<b>Guidance and manufacturer's declaration - electromagnetic emissions</b>		
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal functions. Therefore, the RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	The equipment is not likely to cause any flicker, as the current flow is approx. constant after the start up.	

<b>Guidance and manufacturer's declaration - electromagnetic immunity</b>			
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601-1-2 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV 100 kHz repetition frequency Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV L-N ±2 kV L-PE; N-PE Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption, and voltage variations on power supply input lines IEC 60601-4-11	U <sub>T</sub> =0%, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°)  U <sub>T</sub> =0%, 1 cycle  U <sub>T</sub> =70% 25/30 cycles (at 0°)  U <sub>T</sub> =0%, 250/300 cycles	U <sub>T</sub> =>95%, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°)  U <sub>T</sub> =>95%, 1 cycle  U <sub>T</sub> =70% (30% dipp U <sub>T</sub> ), 25(50Hz)/30(60Hz) cycles (at 0°)  U <sub>T</sub> =>95%, 250(50Hz)/300(60Hz) cycle	Mains power quality should be that of a typical commercial or hospital environment. The device stops and restarts automatically at each dip. In this case does not occur unacceptable pressure drop.
NOTE: U <sub>T</sub> is the A.C. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}, 80 \text{ MHz to } 800 \text{ MHz}$ $d=2,3\sqrt{P}, 800 \text{ MHz to } 2,7 \text{ GHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	<p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p>
Proximity fields from RF wireless communications equipment IEC 61000-4-3	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup>, should be less than the compliance level in each frequency range <sup>b</sup>. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p><sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.</p> <p><sup>b</sup> Over the frequency range 150 kHz to 80 kHz, field strengths should be less than 3 V/m.</p>			

## 4. FUNCTION DESCRIPTION

### 4.1. Function Description

After switching-on the main switch (5) into the position „I“ the network indicator (3) goes on. Taking the suction tubing (6) off the holder (2) will actuate the suction aggregate (9) and the underpressure occurs at the suction mouthpiece (1). After repeated putting the suction tubing into the holder, the suction aggregate turns off. When the separation vessel (11) is filled with waste products, the suction aggregate turns off and the indicator for the separation vessel (4) filling goes on. Then it is necessary to return the suction tubing back into the holder and to empty the separation vessel. During a longer work, mainly with the desalivating mouthpiece, the box temperature may increase; the cooling fan is then automatically actuated. The fan turns off automatically, when the box temperature drops

### 4.2. Detailed Description of the Suction Part Function

The underpressure air flows together with sucked waste products from the oral cavity through the tubing system from the suction mouthpiece (1) at first through the inlet sieve (7), wherein the solid impurities are entrapped. The underpressure air, together with the sucked waste products, free of solid impurities greater than 2 mm, then flows into the separation vessel (11) wherein it is separated from the underpressure air and entrapped into the separation vessel.

Suction air flows through the filter (18) and into the suction unit (9), from which it is pushed through a noise muffler. Therein the output air passes through the output pre-filter (15) and the bacteriologic output filter (14). After passing the filters, the air free of impurities is blown off into a free space under the dental aspirator.

## ASSEMBLY

## 5. STORAGE AND TRANSPORT CONDITIONS

The dental aspirator is delivered ex-work packed in a transport carton, preventing the equipment from being damaged during the transport.

**For the transport, always use an original product packing, if possible.  
Transport the dental aspirator in an upright position.**

Protect the aspirator during transport against moisture, impurity or extreme temperatures. The aspirator in an original packing can be stored in warm, dry and dustfree rooms.

Keep the packing material for possible future return shipments. If this is not possible, dispose the packing material in an environmental-friendly way. The transport carton can be disposed off as the old paper.

**Mobile dental aspirator may only be transported with emptied separation vessel.  
Always empty the separation vessel contents before any transport.**

**Storing or shipping the equipment in any conditions other than those specified below is prohibited.**

### **Climatic conditions for the storage and transport**

**Temperature** –25°C up to +55°C, 24 h up to +70°C

**Relative air humidity** 10% up to 90 % (without condensation)

## 6. INSTALLATION AND INITIAL START-UP

**Before the initial start-up, remove all locking means serving for fixing the equipment during the transport.**

**First actuation can be done only by professional specialist.**

**ANY MODIFICATION OF THIS EQUIPMENT IS FORBIDDEN!**

**This equipment cannot be used nearby other instruments. If this equipment is used nearby other instruments, the equipment must be observed in order to verify normal operations in the configuration it will be used.**

**Instruments may be affected electro-magnetically!**

### 6.1. Environmental Conditions

- The equipment may only be installed and operated in dry, well ventilated and dust-free rooms.
- The mobile dental aspirator should be installed with regard to an easy access for the operation and maintenance and good accessible rating plate.
- The equipment should stand on a plane, sufficiently stable base (be aware of the aspirator weight, see item 3. Technical data).

**Neither the lines for the connection to the electric network, nor the suction tubings may be cranked.**

- To guarantee a trouble-free operation of the aspirator, the room temperature should never drop under +5°C and exceed +40°C. The ideal ambient temperatures are between +10°C up to +25°C.
- Most of the electric energy used by the suction aggregate (9) is converted to heat and given off to the ambient. During a longer lasting work, above all with the desalivating mouthpiece, the box temperature increases over 40°C causing an automatic actuation of the cooling fan. After the room is cooled down under approximately 32°C, the fan turns off again.

### 6.2. Electric Connection

The dental aspirator model 230V is delivered with the safety plug. It is unavoidable to observe the rules of the local electricity works. The network voltage and frequency must be in compliance with the rating plate data.

- In case that the equipment is fix connected to the supply of electric voltage, a disconnecting device must be situated near the equipment.
- If the equipment is connected to the supply of electric voltage by means of a plug, the socket must be good accessible for safety reasons, so that the equipment can safely be unplugged from the network in case of a danger.
- The network distribution must be protected by maximum 10 A.

**The power cord can only be replaced by maintenance personnel!**

### 6.3. Initial Start-up

After unpacking, put the aspirator on the floor, open the door and check whether the separator's cover (10) with sensing electrodes is attached to the separation vessel (11). In case of need, fix the cover to the vessel and then follow the instructions in the Chapter No. 9. Close the door and connect the aspirator's power cord to the socket. Change the main switch (5) over into the position "I", the indicating lamp for the network (3) goes on and the aspirator is ready for the operation.

**Suction cannulas must comply with local legislative and technical requirements for dental suction cannulas and their connection dimensions must be compatible with suction cannulas connectors of 11mm and 16mm diameter**

**Cannulas are not delivered in a sterile state!  
Suction cannulas must undergo complete preparation cycles before first use and after each patient or use following instructions in Chapter 9.7**

## USE

### 7. OPERATION

**The equipment may only be operated by the trained staff!**

**Disconnect from the network in case of a danger (unplug) ! Hot surfaces of the suction aggregate. Do not touch - Danger of burn.**

#### 7.1. Switching on the Dental Aspirator

Switch on the dental aspirator by turning the power-supply switch (5) in the position „I“. The aspirator is put into operation automatically by taking the suction tubing (6) with a mouthpiece (1) off the holder (2). It works uninterruptedly until both tubings are put on the holders again or until the separation vessel (11) is full (during an uniform operation it is filled within about 6 – 10 hours). Suction tips holder is equipped with vacuum regulation of suction unit depends on dentist needs. The state of filling of the separation vessel is evaluated by the separation automatics which interrupts the suction by the aggregate (9) and filling of the separation vessel is indicated automatically with the indicator lamp (4) going on. At this point it is necessary to put the tubings with suction mouthpieces into the holder and switch the power-supply switch off. Afterwards open the door by pulling the grips situated on the side walls, uncouple rubber fixing lugs from the separator cover (10) and pull the separator's vessel out. Hang the separator's cover into the holder on the left side (8). Empty the separation vessel content into the sewer, rinse the vessel with water, connect with the separator's cover in a reversal order. Set the separator into a bowl (12) in the aspirator box (the bowl is detachable). Check the connection of the cover with the vessel, seating of the separation vessel in the bowl and close the door.

#### 7.2. Use of antifoaming tablets

Under specific exhaust conditions, increased exhaust condensate foaming may occur due to the switching off of the exhaust unit when the separation tank is not full. Antifoaming tables (17) must be inserted in order to prevent this from occurring (see the basic accessories); add these tables to the intake screen (7). The tablets will then gradually dissolve in the input screen and will significantly decrease exhaust condensate foaming; the tablets also function as disinfection agents.

#### 7.3. Filter use in the separation vessel

The filter in the separation vessel is an additional layer of protection for the suction unit to prevent moisture from the foam from entering the suction unit. Foam is generated in the separation vessel under specific suction conditions. Foam-blocking tablets are used to eliminate such foam, see Chapter

***Schedule for adding tablets:*** Insert 1 or 2 antifoaming tablets into the intake screen (7) at the end of work and after cleaning and disinfecting the equipment.

**Suction cannulas must undergo complete preparation cycles before first use and after each patient or use following instructions in Chapter 9.7**

**Used suction cannulas intended for repeated use must be cleaned and steam sterilized after every patient (see Chapter 9.7).  
The cannulas are to be exchanged after maximum of 100 cycles of steam sterilizing.  
Cannulas must be visually inspected for any signs of wear or tear before they are used. -Never use the cannulas if it is damaged: discard it immediately.**

It's forbidden to cover all blow-holes placed on upper sides of unit.

Aspirator contains no spare power supply!

Before each switch on running unit it's necessary to check if door of box is closed.

**8. MAINTENANCE INTERVALS - USER / TECHNICIAN**

**Notice!**

The operating entity is obliged to ensure that all tests of the equipment are carried out repeatedly at least once within every 24 months (EN 62353) or in intervals as specified by the applicable national legal regulations. A report must be prepared on the results of the tests (e.g.: according to EN 62353, Annex G), including the measurement methods used.

Maintenance to be done	Chapter	Time interval	To be performed by
Use of antifoaming tablets	7.2	at need	user
Inlet Sieve Cleaning	9.1	per day	user
Lubricating seals and sliding clamps	9.2	Every 15 days	User
Disinfection of Tubings and Separation Vessel	9.3	per day	user
Cleaning and disinfection of the exterior surfaces of the product	9.4	at need	user
Replacement of output pre-filter	9.5	every 3 months	user
Replacement of output filter	9.6	per annum	user
Exchange of suction cannula	9.7	max. 100 steam cleaning cycles	user
Replacing the filter in the separation vessel	9.8	Every 18 months	User
Perform "Repeated Test" according to EN 62353	8	1x in 2 years	qualified technician

**9. MAINTENANCE, CLEANING AND DISINFECTION**

Repair works exceeding the scope of a general maintenance may only be carried out by a qualified specialist or by the manufacturer's customer service. Use only spare parts and accessories permitted by the manufacturer.

Before making any maintenance, repair or cleaning, it is unavoidable to disconnect the equipment from the network (unplug).

Allow the equipment to cool down before maintenance and service!

Reconnect the grounding lead if removed during service work to its previous location once the work is complete.

During the preparation phase for suction cannulas, always wear EC marked personal protective equipment (PPE).

Following works must be carried out from the hygiene and the correct aspirator function points of view.

### 9.1. Inlet Sieve Cleaning

Solid particles (sucked together with a liquid component and the air) are entrapped in the inlet sieve (7) during the equipment operation (Enclosure No. 1, Fig. 1) which therefore must be regularly cleaned, by section 8 (though, always at the end of each shift).

The inlet sieve cleaning is possible only after the moisture inside suction tubings (6) and a sieve is reduced by means of the several seconds lasting air suction through tubings taken off the holder (2). Then turn the power-supply switch (5) off in the position „O“ and slightly lift the cover of the inlet sieve. Afterwards pull the sieve out of its place gripping the holder, remove solid particles and clean the sieve. Then put the sieve on the previous place into the holder and carry out the assembly in the contrary manner to the dismantling.

If the suction of amalgam particles is expected, it is necessary to empty the inlet sieve content into a closing vessel and forward it to the collecting centres.

### 9.2. Lubricating seals and sliding closures

O-rings (sealing rings) and sliding terminal clamps (Enclosure No. 1- pos.1) must be lubricated using a suitable silicone oil for dental applications (e.g. Lubri-Jet spray) every 15 working days.

### 9.3. Disinfection of Tubings and Separation Vessel

For the separation vessel (11) disinfection it is necessary to empty the same, rinse with warm water and clean mechanically with warm water and a disinfectant.

The disinfection of internal suction tubings and of the separation vessel should be performed by section 8, though always at the end of a shift by suction of warm water containing a disinfectant having a limited foaming power through both suction tubings (6) separately. The used disinfectant must be approved in compliance with a valid national legislation on usage for disinfected surface and a type of material. For the disinfectant use, follow the manufacturer's instructions.

### 9.4. Cleaning and disinfection of the exterior surfaces of the product

Use neutral agents for cleaning and disinfection of the external surfaces of the product.

**The use of aggressive cleaning agents and disinfectants comprising alcohol and chlorides may lead to the damage of the surface and the discoloration of the product.**

To find out whether the aspirator works correctly, the following works should be carried out in the defined maintenance intervals:

### 9.5. Replacement of the Output Pre-filter

It is necessary to replace the output pre-filter (15) by section 8. Before making any replacements turn the aspirator power-supply switch (5) off. The pre-filter together with the filter (14) are situated on the left lower part of the aspirator's bogie having the symbol .

For its dismantling, catch the filter holder (13) with one hand and with the other hand slew the fixing lugs (16) by 90° by simultaneous pulling down. Remove the filter holder together with the output filter. Then release the output pre-filter from the cavity wherein also the output filter was situated (check also the state of pollution of the output filter).

For the reassembly, set the output filter into the filter holder, put the output pre-filter onto the filter (direct the output pre-filter with the gummed cloth towards the filter) and insert all the assembly into the cavity space from below. Slew the fixing lugs back by 90° so that they engage into the recessed parts of the filter holder.

**9.6. Replacement of Output Filter**

During the regular use of the device, the output filter (14) must be replaced according to the interval in Chapter 8. The procedure of disassembly and assembly of the filter is the same as for the replacement of output prefilter (15).

**9.7. Cannulas handling and preparation operation**

The preparation instructions have been drafted by supplier of cannulas in accordance with EN ISO 17664 standard requirements.


**During the preparation phase for suction cannulas, always wear EC marked personal protective equipment (PPE).**

TYPE / MANUFACTURER OF CANNULAS DELIVERED FOR THE PRODUCT	Type: Monoart aspirator tips EM21 22910103, Monoart aspirator tips EM21 EVO 22931001 Euronda, Via Chizzalunga 1, 36066 Sandrigo, VI, Italia Tel.(+39) 0444 656185 www.euronda.com
WARNINGS	<ul style="list-style-type: none"> <li>• Store brand new or unused cannulas in a dry, clean and protected environment</li> <li>• Only use EC marked detergents and disinfectants, suitable for use with polypropylene (PP), following the indications concerning concentration, temperature and time of action provided in the manufacturer’s instructions.</li> <li>• Collect and clean the cannulas immediately after use. Dried or sticky residue can make cleaning operations difficult or ineffective.</li> <li>• Never use metal brushes to clean such products.</li> <li>• Damaged cannulas must not be reused.</li> </ul>
TRANSPORTATION AFTER USE	Place the product in a specific container with disinfectant solution according to the concentrations and soak times indicated by the manufacturer.
PREPARATION FOR CLEANING	If necessary, remove visible dirt from the product and then rinse under running water.
MANUAL CLEANING	<ul style="list-style-type: none"> <li>• Use a nylon brush and detergent/disinfectant, suitable for use with polypropylene (PP), to clean the product until all dirt has been removed. Use a bottle brush to clean the insides.</li> <li>• Rinse thoroughly under running water (of good microbiological quality) until all detergent/disinfectant residue has been removed.</li> <li>• Dry well using compressed air.</li> </ul>
ULTRASONIC CLEANING	<ul style="list-style-type: none"> <li>• Place the product in a perforated basket and place it inside the ultrasonic bath filled with detergent or a disinfectant solution suitable for polypropylene (PP).</li> <li>• Set the wash cycle following the manufacturer’s instructions; it is not recommended to set the temperature above 45°C.</li> <li>• Rinse thoroughly under running water (of good microbiological quality) until all detergent residue has been removed.</li> <li>• Check there is no dirt residue on the product; repeat the cleaning cycle whenever necessary.</li> <li>• Dry well using compressed air.</li> </ul> <p>Note: ultrasonic cleaning can be used in conjunction with manual and automatic cleaning and disinfection.</p>
THERMAL DISINFECTOR CLEANING	Use cleaning and disinfectant equipment that is EN ISO 15883 standard-compliant. <ul style="list-style-type: none"> <li>• Position the product using the supplied holders so that the water flows easily and rinses the inner surfaces.</li> </ul>

	<ul style="list-style-type: none"> <li>Follow the manufacturer’s instructions concerning selection of cleaning/disinfection cycles and detergents to be used.</li> <li>On completing the cycle, check there is no dirt residue on the product; repeat the cleaning cycle whenever necessary.</li> </ul>
OPERATING CHECKS AND INSPECTION	Visually inspect the condition and cleanliness level of the product. Discard all damaged or worn products.
PACKAGING	Select EC marked packaging that meets the requirements of Council Directive 93/42/CEE and addenda, suitable for the product and the sterilization procedure, leaving enough space to prevent tension when sealed.
STEAM STERILIZATION	<p>Fractionated pre-vacuum steam sterilization using EN 13060 or EN 285 standard compliant equipment and validated sterilization procedure in accordance with EN ISO 17665-1 standard requirements.</p> <ul style="list-style-type: none"> <li>Use type B cycles.</li> <li>Sterilization temperature: 134°C.</li> <li>Stay time: 5 minutes.</li> <li>Dry time: 10 minutes.</li> </ul>
MAXIMUM DURATION	100 cycles
STORAGE	Storage times shall depend on the type of packaging used and the storage environment. Packaging must be protected against dust, humidity and all contamination risks.

**9.8. Replacing the filter in the separation vessel**

**10. PUTTING OUT OF OPERATION**

When the dental aspirator is not used for a longer time, it is necessary to carry out the cleaning and disinfection of all parts as described in the Articles No. 9.1, 9.3 and then turn the power-supply switch (5) on in the position “I”, take the suction tubings (6) off the holder (2) and let the air intake flow through (for approximately 15-20 min) so that the suction system is fully dried out. Afterwards put the tubings into the holder and turn off the power-supply switch in the position “0”, disconnect the equipment from the network (unplug) and remove the suction mouthpieces (1) from suction tubings.

**DISPOSAL**

**11. DISPOSAL OF APPLIANCE**

- Disconnect the appliance from mains.
- To follow the rules of personal hygiene for works with contaminated material
- Clean device according, section 9.

- To separate, label, packing and providing for decontamination of contaminated parts by course of national regulations.
- Liquidate dental suction unit according to local effective statutes.

**Inside parts of unit can be contaminated with biological material by reason incorrect using. Before clearing and waste disposal pass on special institution for decontamination.**

**TROUBLESHOOTING**

**12. INFORMATION ON REPAIR SERVICE**

Guarantee and post-guarantee repairs are ensured by manufacturer or organizations and repairers denoted by the supplier.

*Warning !*

Manufacturer reserves the right to make changes on the appliance that however will have no impact on substantial characteristics of the appliance.

**13. SOLVING COMMON PROBLEMS**

**Disconnect the appliance from the network before making any interventions therein.**

Works in connection with the removal of failures may only be carried out by a trained stomer service technician.

In case of repair of parts of device which might be contaminated please follow bellow mentioned instruction:

- Disconnect the appliance from mains
- To follow the rules of personal hygiene for works with contaminated material
- Clean device according, section 9
- To separate, label, packing and providing for decontamination of contaminated parts by course of national regulations
- To perform the repair of damaged parts

<b>Failure</b>	<b>Problem and possible reasons</b>	<b>Way of removal</b>
Aspirator does not work	Mains voltage absent	Check mains voltage in the socket, activate a protective element in the electric distribution (fuse, circuit breaker)
	Interrupted supply of el. energy	
	Power-supply switch off	Check the power-supply switch position („I“), switch on the power-supply switch (network indicator goes on)
	Interrupted power cord	Replace the defective part
	Loose clamp of terminal board	Tighten the clamps
	Interrupted motor wiring, defective thermal protection	Replace suction aggregate
	Defective mains fuse of the equipment	Check the state of equipment mains fuse, replace defective fuse (fuses 2 x T 6.3 A are situated in the upper, narrowed part of the aspirator)
	Full separation vessel	Empty the separation vessel content
	High current consumption	Defective starting capacitor, replace the capacitor

	Suction unit is/was awash by liquid	Suction unit pass on repair to service – keep all rools and norms for work with contamination material
Aspirator is switching without taking suction tubings off	Defective control automatics	Replace the defective part
	Defective microswitch in the suction tubings holder	Replace the defective part
Aspirator is noisy	Box door improper closed	Close the box door
	Defective bearings of the suction aggregate	Replace the defective bearing for an identical type
	Incorrectly mounted filter holder	Mount the filter holder into a correct position
	Incorrectly mounted cover of the separation vessel	Mount the separation vessel cover into a correct position
	Suction unit is/was overflowen by liquid.	Suction unit pass on repair to service – keep all rools and norms for work with contamination material
Aspirator efficiency reduced	Heavy polluted output pre-filter	Replace the pre-filter (see Maintenance intervals art.No.8)
	Heavy polluted output filter	Replace the filter (see Maintenance intervals art.No. 8)
	Leaks in the suction system	Check connections, seal untight connections
	Defective control valve in the suction tubings holder	Replace the defective part
	Heavy polluted inlet sieve	Clean the inlet sieve

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
**17. INSTALLATION RECORD**

<b>1. Product:</b> (model)  DO M		<b>2. Serial number:</b>	
<b>3.1. User's name:</b>			
<b>3.2. Address of installation:</b>			
<b>4. Equipment connected to the compressor:</b>			
<b>5. Installation / Commissioning:</b>		<b>6. Contents of operator training:</b>	
A. Product completeness check **	Y	A. Description of the product and functions**	Y
	N		N
B. Documentation completeness check **	Y	B. Product operation: turning on/off, controls, control procedures, data on the display panel, alarms, operation in alarm conditions**	Y
	N		N
C. Installation/connection to equipment **	Y	C. Product maintenance: maintenance intervals, maintenance procedure, service intervals, operating activities**	Y
	N		N
D. Functional test **	Y	D. Safety measures, warnings – their meaning and compliance **	Y
	N		N
<b>Notes::</b>			
<b>7. Operator instructed on safety measures, operations and maintenance:</b>			
Name :		Signature:	
Name:		Signature:	
Name :		Signature:	
<b>8. Installation and instruction performed by –</b> First name/Last name		Signature:	
Company:		Address:	
Phone:			
Email:		Date:	
<b>9. Distributor:</b>			
Company:		Address:	
Contact person:			
Phone:		Email: :	

\*\* mark with an "X" in points 5 and 6 (Y - yes /N - no). Enter any observations from points 5 and 6 into the "Notes" section

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