



DK50 4VR/50

Pre štyri stomatologické súpravy
Do czterech unitów stomatologicznych
Для четырех стоматологических установок



Príslušenstvo/ Akcesoria/ Принадлежности: **BOX MD KJF FS REG AOK**

	DK50 4VR/50	DK50 4VR/50/M	DK50 4VR/50S	DK50 4VR/50S/M
n=1	4	3 - 4	4	3 - 4
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(6-8bar)	250 l.min ⁻¹ / 6bar	200 l.min ⁻¹ / 6bar	250 l.min ⁻¹ / 6bar	200 l.min ⁻¹ / 6bar
(8-10bar)	210 l.min ⁻¹ / 8bar	140 l.min ⁻¹ / 8bar	210 l.min ⁻¹ / 8bar	140 l.min ⁻¹ / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 77 dB(A)	≤ 77 dB(A)	≤ 58 dB(A)	≤ 58 dB(A)

DK50 2x4VR/110/M

Pre rozvody stlačeného vzduchu
Sprężarki do rozprowadzania sprężonego powietrza
Источник для систем распределения воздуха



Príslušenstvo/ Akcesoria/ Принадлежности: **BOX MD KJF FS REG AOK**

	DK50 2x4VR/110/M	DK50 2x4VR/110S/M
n=1	8+	8+
DRYER	✓	✓
BOX	-	✓
(6-8bar)	430 l.min ⁻¹ / 6bar	430 l.min ⁻¹ / 6bar
(8-10bar)	330 l.min ⁻¹ / 8bar	330 l.min ⁻¹ / 8bar
S1-S3	S1 - 100%	S1 - 100%
	≤ 77	≤ 58